

# RNF19A Antibody, Biotin conjugated



PACO56177

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

E3 ubiquitin-protein ligase which accepts ubiquitin from E2 ubiquitin-conjugating enzymes UBE2L3 and UBE2L6 in the form of a thioester and then directly transfers the ubiquitin to targeted substrates, such as SNCAIP or CASR. Specifically ubiquitinates pathogenic SOD1 variants, which leads to their proteasomal degradation and to neuronal protection.

**Gene ID:**

RNF19A

**Uniprot**

Q9NV58

**Synonyms:**

E3 ubiquitin-protein ligase RNF19A (EC 2.3.2.27) (Double ring-finger protein) (Dorfin) (RING finger protein 19A) (p38), RNF19A, RNF19

**Immunogen:**

Recombinant Human E3 ubiquitin-protein ligase RNF19A protein (741-825AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A