

RHCG Antibody, Biotin conjugated



PACO56193

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Functions as an electroneutral and bidirectional ammonium transporter. May regulate transepithelial ammonia secretion.

Gene ID:

RHCG

Uniprot

Q9UBD6

Synonyms:

Ammonium transporter Rh type C (Rh glycoprotein kidney) (Rhesus blood group family type C glycoprotein) (Rh family type C glycoprotein) (Rh type C glycoprotein) (Tumor-related protein DRC2), RHCG, C15orf6 CDRC2 PDRC2 RHGK

Immunogen:

Recombinant Human Ammonium transporter Rh type C protein (416-479AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A