

MMP28 Antibody, Biotin conjugated



PACO56197

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Can degrade casein. Could play a role in tissues homeostasis and repair.

Gene ID:

MMP28

Uniprot

Q9H239

Synonyms:

Matrix metalloproteinase-28 (MMP-28) (EC 3.4.24) (Epilysin), MMP28, MMP25

Immunogen:

Recombinant Human Matrix metalloproteinase-28 protein (351-475AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A