

FASTKD2 Antibody, Biotin conjugated



PACO56201

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Plays an important role in assembly of the mitochondrial large ribosomal subunit.

Gene ID:

FASTKD2

Uniprot

Q9NYY8

Synonyms:

FAST kinase domain-containing protein 2, mitochondrial, FASTKD2, KIAA0971

Immunogen:

Recombinant Human FAST kinase domain-containing protein 2, mitochondrial protein (117-203AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A