

HRH4 Antibody, Biotin conjugated



PACO56221

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

The H4 subclass of histamine receptors could mediate the histamine signals in peripheral tissues. Displays a significant level of constitutive activity (spontaneous activity in the absence of agonist).

Gene ID:

HRH4

Uniprot

Q9H3N8

Synonyms:

Histamine H4 receptor (H4R) (HH4R) (AXOR35) (G-protein coupled receptor 105) (GPRv53) (Pfi-013) (SP9144), HRH4, GPCR105

Immunogen:

Recombinant Human Histamine H4 receptor protein (204-292AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A