

WASF3 Antibody, FITC conjugated



PACO56232

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Downstream effector molecules involved in the transmission of signals from tyrosine kinase receptors and small GTPases to the actin cytoskeleton. Plays a role in the regulation of cell morphology and cytoskeletal organization. Required in the control of cell shape.

Gene ID:

WASF3

Uniprot

Q9UPY6

Synonyms:

Wiskott-Aldrich syndrome protein family member 3 (WASP family protein member 3) (Protein WAVE-3) (Verprolin homology domain-containing protein 3), WASF3, KIAA0900 SCAR3 WAVE3

Immunogen:

Recombinant Human Wiskott-Aldrich syndrome protein family member 3 protein (368-502AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A