

NFE2L3 Antibody, HRP conjugated



PACO56255

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Activates erythroid-specific, globin gene expression.

Gene ID:

NFE2L3

Uniprot

Q9Y4A8

Synonyms:

Nuclear factor erythroid 2-related factor 3 (NF-E2-related factor 3) (NFE2-related factor 3) (Nuclear factor, erythroid derived 2, like 3), NFE2L3, NRF3

Immunogen:

Recombinant Human Nuclear factor erythroid 2-related factor 3 protein (132-260AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A