

# NFE2L3 Antibody, Biotin conjugated



PACO56257

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Activates erythroid-specific, globin gene expression.

**Gene ID:**

NFE2L3

**Uniprot**

Q9Y4A8

**Synonyms:**

Nuclear factor erythroid 2-related factor 3 (NF-E2-related factor 3) (NFE2-related factor 3) (Nuclear factor, erythroid derived 2, like 3), NFE2L3, NRF3

**Immunogen:**

Recombinant Human Nuclear factor erythroid 2-related factor 3 protein (132-260AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A