

AFF2 Antibody, HRP conjugated



PACO56311

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

RNA-binding protein. Might be involved in alternative splicing regulation through an interaction with G-quartet RNA structure.

Gene ID:

AFF2

Uniprot

P51816

Synonyms:

AF4/FMR2 family member 2 (Fragile X E mental retardation syndrome protein) (Fragile X mental retardation 2 protein) (FMR2P) (Protein FMR-2) (Protein Ox19), AFF2, FMR2 OX19

Immunogen:

Recombinant Human AF4/FMR2 family member 2 protein (575-705AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A