

# ARHGEF11 Antibody, HRP conjugated



PACO56315

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May play a role in the regulation of RhoA GTPase by guanine nucleotide-binding alpha-12 (GNA12) and alpha-13 (GNA13). Acts as guanine nucleotide exchange factor (GEF) for RhoA GTPase and may act as GTPase-activating protein (GAP) for GNA12 and GNA13. Involved in neurotrophin-induced neurite outgrowth.

**Gene ID:**

ARHGEF11

**Uniprot**

O15085

**Synonyms:**

Rho guanine nucleotide exchange factor 11 (PDZ-RhoGEF), ARHGEF11, KIAA0380

**Immunogen:**

Recombinant Human Rho guanine nucleotide exchange factor 11 protein (1221-1318AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A