

MSI1 Antibody, Biotin conjugated



PACO56341

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

RNA binding protein that regulates the expression of target mRNAs at the translation level. Regulates expression of the NOTCH1 antagonist NUMB. Binds RNA containing the sequence 5'-GUUAGUUAGUUAGUU-3' and other sequences containing the pattern 5'-[GA]U(1-3)AGU-3'. May play a role in the proliferation and maintenance of stem cells in the central nervous system.

Gene ID:

MSI1

Uniprot

O43347

Synonyms:

RNA-binding protein Musashi homolog 1 (Musashi-1), MSI1

Immunogen:

Recombinant Human RNA-binding protein Musashi homolog 1 protein (221-347AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A