

PEAR1 Antibody, HRP conjugated



PACO56355

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

When overexpressed, reduces the number of both early and late non-adherent myeloid progenitor cells.

Gene ID:

PEAR1

Uniprot

Q5VY43

Synonyms:

Platelet endothelial aggregation receptor 1 (hPEAR1) (Multiple epidermal growth factor-like domains protein 12) (Multiple EGF-like domains protein 12), PEAR1, MEGF12

Immunogen:

Recombinant Human Platelet endothelial aggregation receptor 1 protein (932-1037AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A