

# PLXNB1 Antibody, HRP conjugated



PACO56359

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Receptor for SEMA4D. Plays a role in RHOA activation and subsequent changes of the actin cytoskeleton. Plays a role in axon guidance, invasive growth and cell migration.

**Gene ID:**

PLXNB1

**Uniprot**

O43157

**Synonyms:**

Plexin-B1 (Semaphorin receptor SEP), PLXNB1, KIAA0407 PLXN5 SEP

**Immunogen:**

Recombinant Human Plexin-B1 protein (664-837AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A