

ABCC8 Antibody, Biotin conjugated



PACO56385

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Subunit of the beta-cell ATP-sensitive potassium channel (KATP). Regulator of ATP-sensitive K(+) channels and insulin release.

Gene ID:

ABCC8

Uniprot

Q09428

Synonyms:

ATP-binding cassette sub-family C member 8 (Sulfonylurea receptor 1), ABCC8, HRINS SUR SUR1

Immunogen:

Recombinant Human ATP-binding cassette sub-family C member 8 protein (729-838AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A