

# RNF168 Antibody, FITC conjugated



PACO56420

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:****Gene ID:**

RNF168

**Uniprot**

Q8IYW5

**Synonyms:**

E3 ubiquitin-protein ligase RNF168 (hRNF168) (EC 2.3.2.27) (RING finger protein 168) (RING-type E3 ubiquitin transferase RNF168), RNF168

**Immunogen:**

Recombinant Human E3 ubiquitin-protein ligase RNF168 protein (194-296AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A