

LRRFIP2 Antibody, Biotin conjugated



PACO56449

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May function as activator of the canonical Wnt signaling pathway, in association with DVL3, upstream of CTNNB1/beta-catenin. Positively regulates Toll-like receptor (TLR) signaling in response to agonist probably by competing with the negative FLII regulator for MYD88-binding.

Gene ID:

LRRFIP2

Uniprot

Q9Y608

Synonyms:

Leucine-rich repeat flightless-interacting protein 2 (LRR FLII-interacting protein 2), LRRFIP2

Immunogen:

Recombinant Human Leucine-rich repeat flightless-interacting protein 2 protein (211-342AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A