

# ARL1 Antibody, Biotin conjugated



PACO56727

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

GTP-binding protein that has very low efficiency as allosteric activator of the cholera toxin catalytic subunit, an ADP-ribosyltransferase. Can activate phospholipase D with very low efficiency. Important for normal function of the Golgi apparatus.

**Gene ID:**

ARL1

**Uniprot**

P40616

**Synonyms:**

ADP-ribosylation factor-like protein 1, ARL1

**Immunogen:**

Recombinant Human ADP-ribosylation factor-like protein 1 protein (2-181AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A