

# CCS Antibody, FITC conjugated



PACO56742

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Delivers copper to copper zinc superoxide dismutase (SOD1).

**Gene ID:**

CCS

**Uniprot**

O14618

**Synonyms:**

Copper chaperone for superoxide dismutase (Superoxide dismutase copper chaperone), CCS

**Immunogen:**

Recombinant Human Copper chaperone for superoxide dismutase protein (145-274AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A