

CCS Antibody, Biotin conjugated



PACO56743

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Delivers copper to copper zinc superoxide dismutase (SOD1).

Gene ID:

CCS

Uniprot

O14618

Synonyms:

Copper chaperone for superoxide dismutase (Superoxide dismutase copper chaperone), CCS

Immunogen:

Recombinant Human Copper chaperone for superoxide dismutase protein (145-274AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A