

# CHN2 Antibody, Biotin conjugated



PACO56751

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

GTPase-activating protein for p21-rac. Insufficient expression of beta-2 chimaerin is expected to lead to higher Rac activity and could therefore play a role in the progression from low-grade to high-grade tumors.

**Gene ID:**

CHN2

**Uniprot**

P52757

**Synonyms:**

Beta-chimaerin (Beta-chimerin) (Rho GTPase-activating protein 3), CHN2, ARHGAP3  
BCH

**Immunogen:**

Recombinant Human Beta-chimaerin protein (1-91AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A