

CLN5 Antibody, HRP conjugated



PACO56761

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Plays a role in influencing the retrograde trafficking of lysosomal sorting receptors SORT1 and IGF2R from the endosomes to the trans-Golgi network by controlling the recruitment of retromer complex to the endosomal membrane. Regulates the localization and activation of RAB7A which is required to recruit the retromer complex to the endosomal membrane.

Gene ID:

CLN5

Uniprot

O75503

Synonyms:

Ceroid-lipofuscinosis neuronal protein 5 (Protein CLN5) [Cleaved into: Ceroid-lipofuscinosis neuronal protein 5, secreted form], CLN5

Immunogen:

Recombinant Human Ceroid-lipofuscinosis neuronal protein 5 protein (42-181AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A