

GPR1 Antibody, Biotin conjugated



PACO56811

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Receptor for the inflammation-associated leukocyte chemoattractant chemerin/RARRES2 suggesting a role for this receptor in the regulation of inflammation. Can act as a coreceptor for HIV-1.

Gene ID:

GPR1

Uniprot

P46091

Synonyms:

G-protein coupled receptor 1, GPR1

Immunogen:

Recombinant Human G-protein coupled receptor 1 protein (308-355AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A