

HEPACAM Antibody, HRP conjugated



PACO56821

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in regulating cell motility and cell-matrix interactions. May inhibit cell growth through suppression of cell proliferation.

Gene ID:

HEPACAM

Uniprot

Q14CZ8

Synonyms:

Hepatocyte cell adhesion molecule (Protein hepaCAM), HEPACAM

Immunogen:

Recombinant Human Hepatocyte cell adhesion molecule protein (265-415AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A