

# HOXB4 Antibody, FITC conjugated



PACO56842

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.

**Gene ID:**

HOXB4

**Uniprot**

P17483

**Synonyms:**

Homeobox protein Hox-B4 (Homeobox protein Hox-2.6) (Homeobox protein Hox-2F), HOXB4, HOX2F

**Immunogen:**

Recombinant Human Homeobox protein Hox-B4 protein (88-166AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A