

HOXB4 Antibody, Biotin conjugated



PACO56843

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.

Gene ID:

HOXB4

Uniprot

P17483

Synonyms:

Homeobox protein Hox-B4 (Homeobox protein Hox-2.6) (Homeobox protein Hox-2F), HOXB4, HOX2F

Immunogen:

Recombinant Human Homeobox protein Hox-B4 protein (88-166AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A