

# KIFAP3 Antibody, FITC conjugated



PACO56866

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Involved in tethering the chromosomes to the spindle pole and in chromosome movement. Binds to the tail domain of the KIF3A/KIF3B heterodimer to form a heterotrimeric KIF3 complex and may regulate the membrane binding of this complex.

**Gene ID:**

KIFAP3

**Uniprot**

Q92845

**Synonyms:**

Kinesin-associated protein 3 (KAP-3) (KAP3) (Smg GDS-associated protein), KIFAP3, KIF3AP SMAP

**Immunogen:**

Recombinant Human Kinesin-associated protein 3 protein (660-792AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A