

LRP2 Antibody



PACO56880

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:500-1:1000

Protein Background:

Acts together with cubilin to mediate HDL endocytosis. May participate in regulation of parathyroid-hormone and para-thyroid-hormone-related protein release.

Gene ID:

LRP2

Uniprot

P98164

Synonyms:

Low-density lipoprotein receptor-related protein 2 (LRP-2) (Glycoprotein 330) (gp330) (Megalin), LRP2


Immunogen:

Recombinant Human Low-density lipoprotein receptor-related protein 2 protein (1211-1369AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



IHC image of PACO56880 diluted at 1:500 and staining in paraffin-embedded human kidney tissue performed on a Leica Bond™ system.