

PRCP Antibody, Biotin conjugated



PACO56931

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Cleaves C-terminal amino acid, linked to proline in peptides such as angiotensin II, III and des-Arg9-bradykinin. This cleavage occurs at acid, c pH, but enzymatic activity is retained with some substrates at neutral pH.

Gene ID:

PRCP

Uniprot

P42785

Synonyms:

Lysosomal Pro-X carboxypeptidase (EC 3.4.16.2) (Angiotensinase C) (Lysosomal carboxypeptidase C) (Proline carboxypeptidase) (Prolylcarboxypeptidase) (PRCP), PRCP, PCP

Immunogen:

Recombinant Human Lysosomal Pro-X carboxypeptidase protein (347-443AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A