

# RHOD Antibody, HRP conjugated



PACO56957

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Involved in endosome dynamics. May coordinate membrane transport with the function of the cytoskeleton. Involved in the internalization and trafficking of activated tyrosine kinase receptors such as PDGFRB. Participates in the reorganization of actin cytoskeleton; the function seems to involve WHAMM and includes regulation of filopodia formation and actin filament bundling. Can modulate the effect of DAPK3 in reorganization of actin cytoskeleton and focal adhesion dissolution.

**Gene ID:**

RHOD

**Uniprot**

O00212

**Synonyms:**

Rho-related GTP-binding protein RhoD (Rho-related protein HP1) (RhoHP1), RHOD, ARHD

**Immunogen:**

Recombinant Human Rho-related GTP-binding protein RhoD protein (1-80AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A