

SERPINB4 Antibody, FITC conjugated



PACO56974

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May act as a protease inhibitor to modulate the host immune response against tumor cells.

Gene ID:

SERPINB4

Uniprot

P48594

Synonyms:

Serpin B4 (Leupin) (Peptidase inhibitor 11) (PI-11) (Squamous cell carcinoma antigen 2) (SCCA-2), SERPINB4, PI11 SCCA2

Immunogen:

Recombinant Human Serpin B4 protein (128-208AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A