

# UTY Antibody, HRP conjugated



PACO57009

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Male-specific histone demethylase that catalyzes trimethylated 'Lys-27' (H3K27me3) demethylation in histone H3. Has relatively low lysine demethylase activity.

**Gene ID:**

UTY

**Uniprot**

O14607

**Synonyms:**

Histone demethylase UTY (EC 1.14.11) (Ubiquitously-transcribed TPR protein on the Y chromosome) (Ubiquitously-transcribed Y chromosome tetratricopeptide repeat protein), UTY, KDM6C

**Immunogen:**

Recombinant Human Histone demethylase UTY protein (1-72AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A