

YPEL5 Antibody, HRP conjugated



PACO57013

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

extracellular region, ficolin-1-rich granule lumen, tertiary granule lumen, neutrophil degranulation.

Gene ID:

YPEL5

Uniprot

P62699

Synonyms:

Protein yippee-like 5, YPEL5

Immunogen:

Recombinant Human Protein yippee-like 5 protein (1-121AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A