

# YPEL5 Antibody, Biotin conjugated



PACO57015

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

extracellular region, ficolin-1-rich granule lumen, tertiary granule lumen, neutrophil degranulation.

**Gene ID:**

YPEL5

**Uniprot**

P62699

**Synonyms:**

Protein yippee-like 5, YPEL5

**Immunogen:**

Recombinant Human Protein yippee-like 5 protein (1-121AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A