

FTL3 Antibody, Biotin conjugated



PACO57031

Product Information

Size:

50ug

Reactivity:

Soybean

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

FTL3

Uniprot

E3NYP3

Synonyms:

Immunogen:

Recombinant Soybean FT-like protein (1-176AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A