

AMMECR1 Antibody, Biotin conjugated



PACO57039

Product Information

Size:	Protein Background:
50ug	nucleus.
Reactivity:	Gene ID:
Human	AMMECR1
Source:	Uniprot
Rabbit	Q9Y4X0
Isotype:	Synonyms:
IgG	AMME syndrome candidate gene 1 protein, AMMECR1
Applications:	Immunogen:
ELISA	Recombinant Human AMME syndrome candidate gene 1 protein (146-248AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A