

MYH4 Antibody, Biotin conjugated



PACO57059

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Muscle contraction.

Gene ID:

MYH4

Uniprot

Q9Y623

Synonyms:

Myosin-4 (Myosin heavy chain 2b) (MyHC-2b) (Myosin heavy chain 4) (Myosin heavy chain IIb) (MyHC-IIb) (Myosin heavy chain, skeletal muscle, fetal), MYH4

Immunogen:

Recombinant Human Myosin-4 protein (896-1045AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A