

PACO57072

Product Information

Size:

50ug

Reactivity:

Xenopus laevis

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Transcription factor that acts synergistically with tal1/scl and lmo2 to specify embryonic dorsal mesoderm to a hematopoietic fate.

Gene ID:

gata1-a

Uniprot

P23767

Synonyms:

GATA-binding factor 1-A (Transcription factor xGATA-1A), gata1-a, gata1a

Immunogen:

Recombinant Xenopus laevis GATA-binding factor 1-A protein (1-359AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A