

TLE4 Antibody, FITC conjugated



PACO57078

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Transcriptional corepressor that binds to a number of transcription factors. Inhibits the transcriptional activation mediated by PAX5, and by CTNNB1 and TCF family members in Wnt signaling. The effects of full-length TLE family members may be modulated by association with dominant-negative AES. Essential for the transcriptional repressor activity of SIX3 during retina and lens development and for SIX3 transcriptional auto-repression.

Gene ID:

TLE4

Uniprot

Q04727

Synonyms:

Transducin-like enhancer protein 4 (Grg-4) (Groucho-related protein 4), TLE4, GRG4 KIAA1261

Immunogen:

Recombinant Human Transducin-like enhancer protein 4 protein (284-394AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A