

THOC5 Antibody, FITC conjugated



PACO57090

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	THOC5
Source:	Uniprot
Rabbit	Q13769
Isotype:	Synonyms:
IgG	THO complex subunit 5 homolog (Functional spliceosome-associated protein 79) (fSAP79) (NF2/meningioma region protein pK1.3) (Placental protein 39.2) (PP39.2) (hTREX90), THOC5, C22orf19 KIAA0983
Applications:	Immunogen:
ELISA	Recombinant Human THO complex subunit 5 homolog protein (2-53AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A