

PFKFB4 Antibody, HRP conjugated



PACO57097

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Synthesis and degradation of fructose 2,6-bisphosphate.

Gene ID:

PFKFB4

Uniprot

Q16877

Synonyms:

6-phosphofructo-2-kinase/fructose-2,6-bisphosphatase 4 (6PF-2-K/Fru-2,6-P2ase 4) (PFK/FBPase 4) (6PF-2-K/Fru-2,6-P2ase testis-type isozyme) [Includes: 6-phosphofructo-2-kinase (EC 2.7.1.105); Fructose-2,6-bisphosphatase (EC 3.1.3.46)], PFKFB4

Immunogen:

Recombinant Human 6-phosphofructo-2-kinase/fructose-2,6-bisphosphatase 4 protein (266-368AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A