

# PFKFB4 Antibody, Biotin conjugated



PACO57099

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Synthesis and degradation of fructose 2,6-bisphosphate.

**Gene ID:**

PFKFB4

**Uniprot**

Q16877

**Synonyms:**

6-phosphofructo-2-kinase/fructose-2,6-bisphosphatase 4 (6PF-2-K/Fru-2,6-P2ase 4) (PFK/FBPase 4) (6PF-2-K/Fru-2,6-P2ase testis-type isozyme) [Includes: 6-phosphofructo-2-kinase (EC 2.7.1.105); Fructose-2,6-bisphosphatase (EC 3.1.3.46)], PFKFB4

**Immunogen:**

Recombinant Human 6-phosphofructo-2-kinase/fructose-2,6-bisphosphatase 4 protein (266-368AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A