

# ARSH Antibody, HRP conjugated



PACO57181

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

endoplasmic reticulum lumen, arylsulfatase activity, glycosphingolipid metabolic process, post-translational protein modification.

**Gene ID:**

ARSH

**Uniprot**

Q5FYA8

**Synonyms:**

Arylsulfatase H (ASH) (EC 3.1.6), ARSH

**Immunogen:**

Recombinant Human Arylsulfatase H protein (449-529AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A