

ARSH Antibody, Biotin conjugated



PACO57183

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

endoplasmic reticulum lumen, arylsulfatase activity, glycosphingolipid metabolic process, post-translational protein modification.

Gene ID:

ARSH

Uniprot

Q5FYA8

Synonyms:

Arylsulfatase H (ASH) (EC 3.1.6), ARSH

Immunogen:

Recombinant Human Arylsulfatase H protein (449-529AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A