

PACO57185

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Regulates Wnt proteins sorting and secretion in a feedback regulatory mechanism. This reciprocal interaction plays a key role in the regulation of expression, subcellular location, binding and organelle-specific association of Wnt proteins. Plays also an important role in establishment of the anterior-posterior body axis formation during development.

Gene ID:

WLS

Uniprot

Q5T9L3

Synonyms:

Protein wntless homolog (Integral membrane protein GPR177) (Protein evenness interrupted homolog) (EVI) (Putative NF-kappa-B-activating protein 373), WLS, C1orf139 GPR177

Immunogen:

Recombinant Human Protein wntless homolog protein (123-232AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A