

# SPTBN1 Antibody, HRP conjugated



PACO57209

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Fodrin, which seems to be involved in secretion, interacts with calmodulin in a calcium-dependent manner and is thus candidate for the calcium-dependent movement of the cytoskeleton at the membrane.

**Gene ID:**

SPTBN1

**Uniprot**

Q01082

**Synonyms:**

Spectrin beta chain, non-erythrocytic 1 (Beta-II spectrin) (Fodrin beta chain) (Spectrin, non-erythroid beta chain 1), SPTBN1, SPTB2

**Immunogen:**

Recombinant Human Spectrin &beta; chain, non-erythrocytic 1 protein (2096-2256AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A