

# PDXDC1 Antibody, FITC conjugated



PACO57246

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Golgi apparatus, intracellular membrane-bounded organelle, cadherin binding.

**Gene ID:**

PDXDC1

**Uniprot**

Q6P996

**Synonyms:**

Pyridoxal-dependent decarboxylase domain-containing protein 1 (EC 4.1.1), PDXDC1, KIAA0251

**Immunogen:**

Recombinant Human Pyridoxal-dependent decarboxylase domain-containing protein 1 protein (681-780AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A