

# RBMS3 Antibody, FITC conjugated



PACO57250

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Binds poly(A) and poly(U) oligoribonucleotides.

**Gene ID:**

RBMS3

**Uniprot**

Q6XE24

**Synonyms:**

RNA-binding motif, single-stranded-interacting protein 3, RBMS3

**Immunogen:**

Recombinant Human RNA-binding motif, single-stranded-interacting protein 3 protein (372-437AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A