

NLRP10 Antibody, HRP conjugated



PACO57261

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Inhibits autoprocessing of CASP1, CASP1-dependent IL1B secretion, PYCARD aggregation and PYCARD-mediated apoptosis but not apoptosis induced by FAS or BID. Displays anti-inflammatory activity. Plays a role in adaptive immunity through control of dendritic cell-mediated transport of antigen to the lymph nodes from peripheral sites. Required for immunity against *C. albicans* infection. Involved in the innate immune response by contributing to proinflammatory cytokine release in response to invasive bacterial infection.

Gene ID:

NLRP10

Uniprot

Q86W26

Synonyms:

NACHT, LRR and PYD domains-containing protein 10 (Nucleotide-binding oligomerization domain protein 8), NLRP10, NALP10 NOD8 PYNOD

Immunogen:

Recombinant Human NACHT, LRR and PYD domains-containing protein 10 protein (589-655AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A