

ABCA12 Antibody, FITC conjugated



PACO57270

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Probable transporter involved in lipid homeostasis.

Gene ID:

ABCA12

Uniprot

Q86UK0

Synonyms:

ATP-binding cassette sub-family A member 12 (ATP-binding cassette transporter 12) (ATP-binding cassette 12), ABCA12, ABC12

Immunogen:

Recombinant Human ATP-binding cassette sub-family A member 12 protein (45-158AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A