

CDH7 Antibody, Biotin conjugated



PACO57303

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Cadherins are calcium-dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types.

Gene ID:

CDH7

Uniprot

Q9ULB5

Synonyms:

Cadherin-7, CDH7, CDH7L1

Immunogen:

Recombinant Human Cadherin-7 protein (146-296AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A